


Multiple Clearance Applicators

Multiple Clearance Applicators combine the accuracy of fixed applicators with the versatility of multiple clearance / gap choices in one unit.


Bar Type

- 440C Stainless steel
- 2 clearances



U Shaped

- 440C Stainless steel
- 2 clearances



Ordering Information

Cat. No.	Description
6957	Double Bar 2" – 2&4 mils
5302	Double Bar 2" – 6&10 mils
5303	Double Bar 3" – 2&4 mils
5304	Double Bar 3" – 6&10 mils
5305	Double Bar 4" – 2&4 mils
5306	Double Bar 4" – 6&10 mils
5307	Double Bar 5" – 2&4 mils
5308	Double Bar 5" – 6&10 mils
5309	Double Bar 6" – 2&4 mils
5310	Double Bar 6" – 6&10 mils
5326	U-Bar 2", 2&4 mils
5327	U-Bar 2", 6&10 mils
5328	U-Bar 3", 2&4 mils
5329	U-Bar 3", 6&10 mils
6948	U-Bar 4", 2 & 4 mils
5331	U-Bar 4", 6&10 mils
5332	U-Bar 5", 2&4 mils
5333	U-Bar 5", 6&10 mils
5334	U-Bar 6", 2&4 mils
5335	U-Bar 6", 6&10 mils

Comes complete with:

Applicator bar
Reuseable storage case

Technical Specifications

Gap Clearance μm (mils)	Film Width mm (inch)	Material
50.8 & 101.6 μm (2 & 4 mils)	50.8 mm (2 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	50.8 mm (2 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	76.2 mm (3 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	76.2 mm (3 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	101.6 mm (4 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	101.6 mm (4 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	127.0 mm (5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	127.0 mm (5 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	152.4 mm (6 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	152.4 mm (6 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	38.1 mm (1.5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	38.1 mm (1.5 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	63.5 mm (2.5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	63.5 mm (2.5 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	88.9 mm (3.5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	88.9 mm (3.5 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	114.3 mm (4.5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	114.3 mm (4.5 in)	Stainless Steel
50.8 & 101.6 μm (2 & 4 mils)	139.7 mm (5.5 in)	Stainless Steel
152.4 & 254.0 μm (6 & 10 mils)	139.7 mm (5.5 in)	Stainless Steel